



Recombinant Human ATP-binding cassette sub-family F member 3 (ABCF3)

Product Code	CSB-YP889115HU
Abbreviation	ABCF3
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9NUQ8
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	ATCAEILRS EFPEIDGQVF DYVTGVLHSG SADFESVDDL VEAVGELLQE VSGDSKDDAG IRAVCQRMYN TLRLAEPQSQ GNSQVLLDAP IQLSKITENY DCGTKLPGLL KREQSSTVNA KKLEKAEARL KAKQEKRSEK DTLKTSNPLV LEEASASQAG SRKESRLESS GKNKSYDVRI ENFDVSFGDR VLLAGADVNL AWGRRYGLVG RNGLGKTTLL KMLATRSRLV PAHISLLHVE QEVAGDDTPA LQSVLESDSV REDLLRRERE LTAQIAAGRA EGSEAAELAE IYAKLEEIEA DKAPARASVI LAGLGFTPCKM QQQPTREFSG GWRMRLALAR ALFARPDLLL LDEPTNMLDV RAILWLENYL QTWPSTILVV SHDRNFLNAI ATDIIHLHSQ RLDGYRGDFE TFIKSKQERL LNQQREYEAQ QQYRQHIQVF IDRFYRNANR ASQVQSKLKM LEKLPKPV DKESEVVMKF PDGFEKFSPP ILQLDEVDFY YDPKHVIFSR LSVSADLESR ICVVGENGAG KSTMLKLLL DLAPVRGIRH AHRNLKIGYF SQHHVEQLDL NVSAVELLAR KFPGRPEEEY RHQLGRYGIS GELAMRPLAS LSGGQKSRVA FAQMTMPCPN FYILDEPTNH LDMETIEALG RALNNFRGGV ILVSHDERFI RLVCRELWVC EGGGVTRVEG GFDQYRALLQ EQFRREGFL
Source	Yeast
Target Names	ABCF3
Protein Names	Recommended name: ATP-binding cassette sub-family F member 3
Expression Region	2-709
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final



concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.