



Recombinant Human Armadillo repeat-containing protein 1 (ARMC1)

Product Code	CSB-BP889123HU
Abbreviation	ARMC1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9NVT9
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MNSSTSTMSE EPDALSVVNQ LRDLAADPLN RRAIVQDQGC LPGLILFMDH PNPPVVHSAL LALRYLAECR ANREKMKGEL GMMLSLQNVI QKTTTPGETK LLASEIYDIL QSSNMADGDS FNEMNSRRRK AQFFLGTTNK RAKTVVLHID GLDDTSRRNL CEEALLKIKG VISFTFQMAV QRCVVRIRSD LKAEALASAI ASTKVMKAQQ VVKSESGEEM LVPFQDTPVE VEQNTELPDY LPEDESPTKE QDKAVSRVGS HPEGGASWLS TAANFLSRSF YW
Source	Baculovirus
Target Names	ARMC1
Protein Names	Recommended name: Armadillo repeat-containing protein 1
Expression Region	1-282
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.