



Recombinant Human Notchless protein homolog 1 (NLE1)

Product Code	CSB-BP889125HU
Abbreviation	NLE1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9NVX2
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	AAAVPDEAV ARDVQRLLVQ FQDEGGQLLG SPFDVPVDIT PDRLQLVCNA LLAQEDPLPL AFFVHDAEIV SSLGKTLESQ AVETEKVLDI IYQPQAIFRV RAVTRCTSSL EGHSEAVISV AFSPTGKYLA SGSGDTTVRF WDLSTETPHF TCKGHRHWVL SISWSPDGRK LASGCKNGQI LLWDPSTGKQ VGRTLAGHSK WITGLSWEPL HANPECRYVA SSSKDGSVRI WDTTAGRCER ILTGHTQSVT CLRWGGDGLL YSASQDRTIK VWRAHDGVLC RTLQGHGHWW NTMALSTDYA LRTGAFEPAE ASVNPQDLQG SLQELKERAL SRYNLVRGQG PERLVSGSDD FTLFLWSPA E DKKPLTRMTG HQALINQVLF SPDSRIVASA SFDKSIKLWD GRTGKYLASL RGHVAAVYQI AWSADSRLLV SGSSDSTLKV WDVKAQKLAM DLPGHAEVY AVDWSPDGQR VASGGKDKCL RIWRR
Source	Baculovirus
Target Names	NLE1
Protein Names	Recommended name: Notchless protein homolog 1
Expression Region	2-485
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.