



# Recombinant Human Pre-mRNA-splicing factor RBM22 (RBM22)

<b>Product Code</b>	CSB-EP889126HU
<b>Abbreviation</b>	RBM22
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q9NW64
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	ATSLGSNTY NRQNWEDADF PILCQTCLGE NPYIRMTKEK YGKECKICAR PFTVFRWCPG VRMRFKKTEV CQTC SKLKNV CQTCLLDLEY GLPIQVRDAG LSFKDDMPKS DVNKEYYTQN MEREISNSDG TRPVGMLGKA TSTSDMLLKL ARTTPYYKRN RPHICSFVVK GECKRGE ECP YRHEKPTDPD DPLADQNIKD RYYGINDPVA DKLLKRASTM PRLDPPEDKT ITTLYVGGLG DTITETDLRN HFYQFGEIRT ITVVQRQQCA FIQFATRQAA EVAAEKSFNK LIVNGRRLNV KWGRSQAARG KEKEKDGTTD SGIKLEPVPVPG LPGALPPPPA AEEEASANYF NLPPSGPPAV VNIALPPPPG IAPPPPPGFG PHMFHPMGPP PPFMRAPGPI HYPSQDPQRM GAHAGKHSSP
<b>Source</b>	E.coli
<b>Target Names</b>	RBM22
<b>Protein Names</b>	Recommended name: Pre-mRNA-splicing factor RBM22 Alternative name(s): RNA-binding motif protein 22 Zinc finger CCCH domain-containing protein 16
<b>Expression Region</b>	2-420
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Target Details</b>	This gene encodes an RNA binding protein. The encoded protein may play a role in cell division and may be involved in pre-mRNA splicing. Related pseudogenes exist on chromosomes 6, 7, 9, 13, 16, 18, and X.
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



## Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.