



# Recombinant Human Uridine-cytidine kinase-like 1 (UCKL1)

<b>Product Code</b>	CSB-MP889133HU
<b>Abbreviation</b>	UCKL1
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q9NWZ5
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MAAPPARADA DPSPTSPTA RDTPGRQAEK SETACEDRSN AESLDRLLPP VGTGRSPRKR TTSQCKSEPP LLRTSKRTIY TAGRPPWYNE HGTQSKEAFA IGLGGGSASG KTTVARMIE ALDVPWVLL SMDSFYKVLV EQQQEQAAHN NFNFDHPDAF DFDLIISTLK KLKQGKSVKV PIYDFTTHSR KKDWKTYLGA NVIIFEGIMA FADKTLELL DMKIFVDTDS DIRLVRLRR DISERGRDIE GVIKQYNKFV KPSFDQYIQP TMLRADIVVP RSGSNTVAID LIVQHVHSQL EERELSVRAA LASAHQCHPL PRTLVLKST PQVRGMHTII RDKETSDEF IFYSKRLMRL LIEHALSFLP FQDCVVQTPQ GQDYAGKCYA GKQITGVSIL RAGETMEPAL RAVCKDVRIG TILIQTNQLT GEPELHYLRL PKDISDDHVI LMDCTVSTGA AAMMAVRVLL DHDVPEDKIF LLSLLMAEMG VHSVAYAFPR VRIITAVDK RVNDLFRIIP GIGNFGDRYF GTDAVPDGSD EEEVAYTG
<b>Source</b>	Mammalian cell
<b>Target Names</b>	UCKL1
<b>Protein Names</b>	Recommended name: Uridine-cytidine kinase-like 1 EC= 2.7.1.48
<b>Expression Region</b>	1-548
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.