



Recombinant Human Uridine-cytidine kinase-like 1 (UCKL1)

Product Code	CSB-YP889133HU
Abbreviation	UCKL1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9NWZ5
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MAAPPARADA DPSPTSPTA RDTPGRQAEK SETACEDRSN AESLDRLLPP VGTGRSPRKR TTSQCKSEPP LLRTSKRTIY TAGRPPWYNE HGTQSKEAFA IGLGGGSASG KTTVARMIE ALDVPWVLL SMDSFYKVLV EQQQEQAAHN NFNFDHPDAF DFDLIISTLK KLKQGKSVKV PIYDFTTHSR KKDWKTYGA NVIIFEGIMA FADKTLELL DMKIFVDTDS DIRLVRRLRR DISERGRDIE GVIKQYNKFV KPSFDQYIQP TMLRADIVVP RGSNTVAID LIVQHVHSQL EERELSVRAA LASAHQCHPL PRTLVLKST PQVRGMHTII RDKETSDEF IFYSKRLMRL LIEHALSFLP FQDCVVQTPQ GQDYAGKCYA GKQITGVSIL RAGETMEPAL RAVCKDVRIG TILIQTNQLT GEPELHYLRL PKDISDDHVI LMDCTVSTGA AAMMAVRVLL DHDVPEDKIF LLSLLMAEMG VHSVAYAFPR VRIITAVDK RVNDLFRIIP GIGNFGDRYF GTDAVPDGSD EEEVAYTG
Source	Yeast
Target Names	UCKL1
Protein Names	Recommended name: Uridine-cytidine kinase-like 1 EC= 2.7.1.48
Expression Region	1-548
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.