



Recombinant Human Actin-related protein 10 (ACTR10)

Product Code	CSB-EP889156HU-B
Abbreviation	ACTR10
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9NZ32
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MPLYEGLGSG GEKTAVVIDL GEAFKCGFA GETGPRCIIP SVIKRAGMPK PVRVVQYNIN TEELYSYLKE FIHILYFRHL LVNPRDRRVV IIESVLCPSH FRETLTRVLF KYFEVPSVLL APSHLMALLT LGINSAMVLD CGYRESLVLP IYEGIPVLNC WGALPLGGKA LHKELETQLL EQCTVDTSVA KEQSLPSVMG SVPEGVLEDI KARTCFVSDL KRGLKIQAAK FNIDGNNERP SPPPNVDYPL DGEKILHILG SIRDSVVEIL FEQDNEEQSV ATLILDSLQI CPIDTRKQLA ENLVVIGGTS MLPGFLHRL AEIRYLVEKP KYKKALGTKT FRIHTPPAKA NCVAWLGGAI FGALQDILGS RSVSKEYYNQ TGRIPDWCSL NNPPLEMMFD VGKTQPPLMK RAFSTEK
Source	E.coli
Target Names	ACTR10
Protein Names	Recommended name: Actin-related protein 10 Alternative name(s): Actin-related protein 11 Short name= hARP11
Expression Region	1-417
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.