



Recombinant Human Interleukin-36 beta (IL36B), partial

Product Code	CSB-EP889159HU-B
Abbreviation	IL36B
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9NZH7
Storage Buffer	Lyophilized from Tris/PBS-based buffer, 6% Trehalose, pH 8.0
Product Type	Recombinant Proteins
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	REAAPKSYAIRDSRQMVWVLSGNSLIAAPLSRSIKPVTLHLIACRDTEFSDKEK GNMVYLGIGKGDLCFLCAEIQGKPTLQLKLQGSQDNIGKDTCWKLVGIHTCINL DVRESCFMGTLDQWGIGVGRKKWKSSFQHHHLRKKDKDFSSMRTN
Source	E.coli
Target Names	IL36B
Protein Names	Recommended name: Interleukin-36 beta Alternative name(s): FIL1 eta Interleukin-1 eta Short name= IL-1 eta Interleukin-1 family member 8 Short name= IL-1F8 Interleukin-1 homolog 2 Short name= IL-1H2
Expression Region	5-157aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Partial
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.