



Recombinant Human Denticleless protein homolog (DTL)

Product Code	CSB-BP889160HU
Abbreviation	DTL
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9NZJ0
Storage Buffer	Lyophilized from Tris/PBS-based buffer, 6% Trehalose, pH 8.0
Product Type	Recombinant Proteins
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MLFNSVLRQP QLGVLRNGWS SQYPLQSLLT GYQCSGNDEH TSYGETGVPV PPFGCTFSSA PNMEHVLAVA NEEGFVRLYN TESQSFRKKC FKEWMAHWNA VFDLAWVPGE LKLVTAAGDQ TAKFWDVKAG ELIGTCKGHQ CSLKSVAFSK FEKAVFCTGG RDGNIMVWDT RCNKKDGFYR QVNQISGAHN TSDKQTPSKP KKKQNSKGLA PSVDFQQSVT VVLFQDENTL VSAGAVDGII KVDWLRKNYT AYRQEPIASK SFLYPGSSTR KLGYSLLILD STGSTLFANC TDDNIYMFNM TGLKTPVAI FNGHQNSTFY VKSSLSPDDQ FLVSGSSDEA AYIWKVSTPW QPPTVLLGHS QEVTSVCWCP SDFTKIATCS DDNTLKIWRL NRGLEEKPGG DKLSTVGWAS QKKKESRPGL VTVTSSQSTP AKAPRAKCNP SNSSPSSAAC APSCAGDLPL PSNTPTFSIK TSPAKARSP NRRGSVSSVS PKPPSSFKMS IRNWVTRTPS SSPPITPPAS ETKIMSPRKA LIPVSQKSSQ AEACSESRNR VKRRLDSSCL ESVKQKCVKS CNCVTELDGQ VENLHLDLCC LAGNQEDLSK DSLGPTKSSK IEGAGTSISE PPSPISPYAS ESCGTLPLPL RPCGEGSEMV GKENSSENK NWLLAMAAKR KAENPSRSP SSQTPNSRRQ SGKKLPSPVT ITPSSMRKIC TYFHRKSQED FCGPEHSTEL
Source	Baculovirus
Target Names	DTL
Protein Names	Recommended name: Denticleless protein homolog Alternative name(s): DDB1- and CUL4-associated factor 2 Lethal(2) denticleless protein homolog Retinoic acid-regulated nuclear matrix-associated protein
Expression Region	1-730
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein

**Reconstitution**

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.