



Recombinant Human Denticleless protein homolog (DTL)

| | |
|--------------------------|---|
| Product Code | CSB-MP889160HU |
| Abbreviation | DTL |
| Storage | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |
| Uniprot No. | Q9NZJ0 |
| Storage Buffer | Lyophilized from Tris/PBS-based buffer, 6% Trehalose, pH 8.0 |
| Product Type | Recombinant Proteins |
| Immunogen Species | Homo sapiens (Human) |
| Purity | ≥85% (SDS-PAGE) |
| Sequence | MLFNSVLRQP QLGVLRNGWS SQYPLQSLLT GYQCSGNDEH TSYGETGVPV PPFGCTFSSA PNMEHVLAVA NEEGFVRLYN TESQSFRKKC FKEWMAHWNA VFDLAWVPG E LKLVTAAAGDQ TAKFWDVKAG ELIGTCKGHQ CSLKSVAFSK FEKAVFCTGG RDGNIMVWDT RCNKKDGFYR QVNQISGAHN TSDKQTPSKP KKKQNSKGLA PSVDFQQSVT VVLFQDENTL VSAGAVDGII KVDWLRKNYT AYRQEPIASK SFLYPGSSTR KLGYSLLILD STGSTLFANC TDDNIYMFNM TGLKTPVAI FNGHQNSTFY VKSSLSPDDQ FLVSGSSDEA AYIWKVSTPW QPPTVLLGHS QEVTSVCWCP SDFTKIATCS DDNTLKIWRL NRGLEEKPGG DKLSTVGWAS QKKKESRPGL VTVTSSQSTP AKAPRAKCNP SNSSPSSAAC APSCAGDLPL PSNTPTFSIK TSPAKARSP NRRGSVSSVS PKPPSSFKMS IRNWVTRTPS SSPPITPPAS ETKIMSPRKA LIPVSQKSSQ AEACSESRNR VKRRLDSSCL ESVKQKCVKS CNCVTELDGQ VENLHLDLCC LAGNQEDLSK DSLGPTKSSK IEGAGTSISE PPSPISPYAS ESCGTLPLPL RPCGEGSEMV GKENSSENK NWLLAMAAR KAENPSRSP SSQTPNSRRQ SGKKLPSPVT ITPSSMRKIC TYFHRKSQED FCGPEHSTEL |
| Source | Mammalian cell |
| Target Names | DTL |
| Protein Names | Recommended name: Denticleless protein homolog Alternative name(s): DDB1- and CUL4-associated factor 2 Lethal(2) denticleless protein homolog Retinoic acid-regulated nuclear matrix-associated protein |
| Expression Region | 1-730 |
| Notes | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week. |
| Tag Info | Tag type will be determined during the manufacturing process. |
| Protein Length | full length protein |

**Reconstitution**

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.