



Recombinant Human MAGUK p55 subfamily member 6 (MPP6)

Product Code	CSB-EP889164HU
Abbreviation	MPP6
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9NZW5
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	<p> MQQVLENLTE LPSSTGAEI DLIFLKGIME NPIVKSLAKA HERLEDSKLE AVSDNNLELV NEILEDITPL INVDEVAEL VGILKEPHFQ SLLEAHDIVA SKCYDSPSS PEMNSSINN QLLPVDAIRI LGIHKRAGEP LGVTFRVENN DLVIARILHG GMIDRQGLLH VGDIKEVNG HEVGNNPKEL QELLKNISGS VTLKILPSYR DTITPQQVFV KCHFDPNPYN DNLIPECKEAG LKFSKGEILQ IVNREDPNWW QASHVKEGGS AGLIPSQFLE EKRKAFVRRD WDNSGPFPGT ISSKKKKMM YLTTRNAEFD RHEIQIYEEV AKMPPFQRKT LVLIGAQGVG RRSLKNRFIV LNPTFRGTTV PFTSRKPRED EKDGQAYKFV SRSEMEADIK AGKYLEHGEY EGNLYGTKID SILEVVQTGR TCILDVNPQA LKVLRTSEFM PYVVFIAAPE LETLRAMHKA VVDAGITTKL LTDSLKKTIV DESARIQRAY NHYFDLIIN DNLDKAFEKL QTAIEKLRME PQWVPISWVY </p>
Source	E.coli
Target Names	MPP6
Protein Names	Recommended name: MAGUK p55 subfamily member 6 Alternative name(s): Veli-associated MAGUK 1 Short name= VAM-1
Expression Region	1-540
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.