



Recombinant Human NADH dehydrogenase [ubiquinone] 1 alpha subcomplex assembly factor 4 (NDUFAF4)

Product Code	CSB-EP889166HU-B
Abbreviation	NDUFAF4
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9P032
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MGALVIRGIR NFNLENRAER EISKMKPSVA PRHPSTNSLL REQISLYPEV KGEIARKDEK LLSFLKDYYV DSKDPVSSLQ VKAAETCQEP KEFRLPKDHH FDMINIKSIP KKGISIVEAL TLLNNHKLFP ETWTAEKIMQ EYQLEQKDVN SLLKYFVTFE VEIFPPEDKK AIRSK
Source	E.coli
Target Names	NDUFAF4
Protein Names	Recommended name: NADH dehydrogenase [ubiquinone] 1 alpha subcomplex assembly factor 4 Alternative name(s): Hormone-regulated proliferation-associated protein of 20 kDa
Expression Region	1-175
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.