



Recombinant Human Forkhead box protein J2 (FOXJ2)

Product Code	CSB-MP889168HU
Abbreviation	FOXJ2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9P0K8
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	ASDLESSLT SIDWLPQLTL RATIEKLGSA SQAGPPGSSR KCSPGSPTDP NATLSKDEAA VHQDGKPRYS YATLITYAIN SSPAKKMTLS EIYRWICDNF PYYKNAGIGW KNSIRHNLSL NKCFRKVPRP RDDPGKGSYW TIDTCPDISR KRRHPPDDDL SQDSPEQEAS KSPRGGVAGS GEASLPPEGN PQMSLQSPTS IASYSQGTGS VDGGAVAAGA SGRESAEGPP PLYNTNHDFK FSYSEINFQD LSWSFRNLYK SMLEKSSSSS QHGFSSLLGD IPPSNYYMY QQQQPPPPQQ QQQQQQPPQP PPQQSQPQQQ QAPAQGPSAV GGAPPLHTPS TDGCTPPGGK QAGAEGYGPP PVMAMHPPPL QHGGYHPHQH HPHSHPAQQP PPPQPAQQQ APINNTGFAF PSDWCSNIDS LKESFKMVNR LNWSSIEQSQ FSELMESLRQ AEQKNWTLDQ HHIANLCDSL NHFLTQTGHV PPQGGTHRPP APARIADSCA LTSGKQESAM SQVNSYGHPQ APHLYPGPSP MYPIPTQDSA GYNRPAHHMV PRPSVPPPGA NEEIPDDFDW DLIT
Source	Mammalian cell
Target Names	FOXJ2
Protein Names	Recommended name: Forkhead box protein J2 Alternative name(s): Fork head homologous X
Expression Region	2-574
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.