



# Recombinant Human Forkhead box protein J2 (FOXJ2)

<b>Product Code</b>	CSB-YP889168HU
<b>Abbreviation</b>	FOXJ2
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q9P0K8
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	ASDLESSLT SIDWLPQLTL RATIEKLGSA SQAGPPGSSR KCSPGSPTDP NATLSKDEAA VHQDGKPRYS YATLITYAIN SSPAKKMTLS EIYRWICDNF PYYKNAGIGW KNSIRHNLSL NKCFRKVPRP RDDPGKGSYW TIDTCPDISR KRRHPPDDDL SQDSPEQEAS KSPRGGVAGS GEASLPPEGN PQMSLQSPTS IASYSQGTGS VDGGAVAAGA SGRESAEGPP PLYNTNHDFK FSYSEINFQD LSWSFRNLYK SMLEKSSSSS QHGFSSLLGD IPPSNYYMY QQQQPPPPQQ QQQQQPPQP PPQQSQPQQQ QAPAQGPSAV GGAPPLHTPS TDGCTPPGGK QAGAEGYGPP PVMAMHPPPL QHGGYHPHQH HPHSHPAQQP PPPQPQAQQQ APINNTGFAF PSDWCSNIDS LKESFKMVNR LNWSSIEQSQ FSELMESLRQ AEQKNWTLQ HHIANLCDSL NHFLTQTGHV PPQGGTHRPP APARIADSCA LTSGKQESAM SQVNSYGHPQ APHLYPGPSP MYPIPTQDSA GYNRPAHHMV PRPSVPPPGA NEEIPDDFDW DLIT
<b>Source</b>	Yeast
<b>Target Names</b>	FOXJ2
<b>Protein Names</b>	Recommended name: Forkhead box protein J2 Alternative name(s): Fork head homologous X
<b>Expression Region</b>	2-574
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



## Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.