



Recombinant Human Cleavage and polyadenylation specificity factor subunit 2 (CPSF2)

Product Code	CSB-YP889185HU
Abbreviation	CPSF2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9P2I0
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	<p> MTSIIKLTTL SGVQEEALC YLLQVDEFRL LLDCGWDEHF SMDIIDSLRK HVHQIDAVLL SHPDPLHLGA LPYAVGKGLL NCAIYATIPV YKMGQMFMYD LYQSRHNTED FTLFTLDDVD AAFDKIQQLK FSQIVNLKGGK GHGLSITPLP AGHMIGGTIW KIVKDGEELI VYAVDFNHRK EIHLNGCSLE MLSRPSSLIT DSNATYVQP RRRQRDEQLL TNVLETLRGD GNVLIAVDTA GRVLELAQLL DQIWRTKDAG LGVYSLALLN NVSYNVVEFS KSQVEWMSDK LMRCFEDKRN NPFQFRHLSL CHGLSDLARV PSPKVVLASQ PDLECGFSRD LFIQWCQDPK NSIILTYRTT PGTARFLID NPSEKITEIE LRKRVKLEGG ELEEYLEKEK LKKEAAKKLE QSKEADIDSS DESDIEEDID QPSAHKTKHD LMMKGEGRK GSFFKQAKKS YPMFPAPEER IKWDEYGEII KPEDFLVPEL QATEEEKSKL ESGLTNGDEP MDQDLSDVPT KCISTTESIE IKARVTYIDY EGRSDGDSIK KIINQMKPRQ LIIVHGPEEA SQDLAECRA FGGKDIKVYM PKLHETVDAT SETHIQVRL KDSLVSSLQF CKAKDAELAW IDGVLDMRVS KVDTGVILEE GELKDDGDS EMQVEAPSDS SVIAQQKAMK SLFGDDEKET GEESEIPTL EPLPPHEVPG HQSVFMNEPR LSDFKQVLLR EGIQAEFVGG VLVCCNNQVAV RRTETGRIGL EGCLCQDFYR IRDLLYEQYA IV </p>
Source	Yeast
Target Names	CPSF2
Protein Names	Recommended name: Cleavage and polyadenylation specificity factor subunit 2 Alternative name(s): Cleavage and polyadenylation specificity factor 100 kDa subunit Short name= CPSF 100 kDa subunit
Expression Region	1-782
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.

**Protein Length**

full length protein

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.