



Recombinant Human Protein IMPACT (IMPACT)

Product Code	CSB-YP889191HU
Abbreviation	IMPACT
Storage	<p>The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.</p> <p>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.</p>
Uniprot No.	Q9P2X3
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MAEGDAGSDQ RQNEEIEAMA AIYGEEWCVI DDCAKIFCIR ISDDIDDPKW TLCLQVMLPN EYPGTAPPIY QLNAPWLKGQ ERADLSNSLE EIIYQNIAGES ILYLWVEKIR DVLIQKSQMT EPGPDVKKKT EEEDVECEDD LILACQPSS LKALDFDISE TRTEVEVEEL PPIDHGIPIT DRRSTFQ AHL APVVCPKQVK MVLSKLYENK KIASATHNIY AYRIYCEDKQ TFLQDCEDDG ETAAGGRL LH LMEILNVKNV MVVVSRWYGG ILLGPDRFKH INNCARNILV EKNYTNSPEE SSKALGKNKK VRKDKKRNEH
Source	Yeast
Target Names	IMPACT
Protein Names	Recommended name: Protein IMPACT Alternative name(s): Imprinted and ancient gene protein homolog
Expression Region	1-320
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	<p>The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.</p> <p>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.</p>