



Recombinant *Schizosaccharomyces pombe* Uncharacterized tRNA-binding protein C30C2.04 (SPAC30C2.04)

Product Code	CSB-EP889215SXV
Abbreviation	SPAC30C2.04
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9P6K7
Product Type	Recombinant Protein
Immunogen Species	<i>Schizosaccharomyces pombe</i> (strain 972 / ATCC 24843) (Fission yeast)
Purity	>85% (SDS-PAGE)
Sequence	MSFLWAKEFC VLKYPYKLF THKFSYLLRF ETVTLNNIRS PQFYAKLTFS HHQPTVSGTM STELKFISKY LQISIPETKE GPVSSLFKAV SEQKPELLGN TDFEKAQILE WTTKAFSPIE TQSIVEQLDE FLKSSTFIAQ DSGISVADLA VYARIHSYIC GLSAKEGYKL NNVCRWDFDI QHQESVMEAA NSMSMKLANI DLNAPKIQRP SVIKKDKKEK KEGKPSQEAS VKSVEKAPKG LEGAKKEKQN KKEKKDKKDK KDKKEKAPKE PPKAATPVPS MIDFRIGFIE KAVKHPNADS LYVSTIHCGD AEGPRTVCSG LVKYIPLEQM QQRKVIVVAN LKPVMNRSVK SQAMVFCASS PDKSVVEFVL PPENAEIGDR LTFEGFDTEE PEAQLNPKRK IWEAIQPGFT SGEDLICGYK DESGLHRLFV KGKKDLGFCK AQTVVNGTLS
Source	E.coli
Target Names	SPAC30C2.04
Protein Names	Recommended name: Uncharacterized tRNA-binding protein C30C2.04
Expression Region	1-450
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.