



Recombinant *Schizosaccharomyces pombe* Uncharacterized tRNA-binding protein C30C2.04 (SPAC30C2.04)

| | |
|--------------------------|---|
| Product Code | CSB-MP889215SXV |
| Abbreviation | SPAC30C2.04 |
| Storage | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |
| Uniprot No. | Q9P6K7 |
| Product Type | Recombinant Protein |
| Immunogen Species | <i>Schizosaccharomyces pombe</i> (strain 972 / ATCC 24843) (Fission yeast) |
| Purity | >85% (SDS-PAGE) |
| Sequence | MSFLWAKEFC VLKYPYKLF THKFSYLLRF ETVTLNNIRS PQFYAKLTFS HHQPTVSGTM STELKFISKY LQISIPETKE GPVSSLFKAV SEQKPELLGN TDFEKAQILE WTTKAFSPIE TQSIVEQLDE FLKSSTFIAQ DSGISVADLA VYARIHSYIC GLSAKEGYKL NNVCRWDFDI QHQESVMEAA NSMSMKLANI DLNAPKIQRP SVIKKDKKEK KEGKPSQEAS VKSVEKAPKG LEGAKKEKQN KKEKKDKKDK KDKKEKAPKE PPKAATPVPS MIDFRIGFIE KAVKHPNADS LYVSTIHCGD AEGPRTVCSG LVKYIPLEQM QQRKVIVVAN LKPVNMRSVK SQAMVFCASS PDKSVVEFVL PPENAEIGDR LTFEGFDTEE PEAQLNPKRK IWEAIQPGFT SGEDLICGYK DESGLHRLFV KGKKDLGFCK AQTVVNGTLS |
| Source | Mammalian cell |
| Target Names | SPAC30C2.04 |
| Protein Names | Recommended name: Uncharacterized tRNA-binding protein C30C2.04 |
| Expression Region | 1-450 |
| Notes | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week. |
| Tag Info | Tag type will be determined during the manufacturing process. |
| Protein Length | full length protein |
| Reconstitution | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference. |
| Shelf Life | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life |



of lyophilized form is 12 months at -20°C/-80°C.