



Recombinant Schizosaccharomyces pombe ATP-dependent rRNA helicase rrp3 (rrp3)

Product Code	CSB-EP889218SXV-B
Abbreviation	rrp3
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9P6N8
Product Type	Recombinant Protein
Immunogen Species	Schizosaccharomyces pombe (strain 972 / ATCC 24843) (Fission yeast)
Purity	>85% (SDS-PAGE)
Sequence	MAPSEKKLTE DKKNSSLNKK IETSNSSEK SSENNGDSQ NNEAPKTFKE LGVIDELCEA CEKLGFKTPT PIQQEAIPVV LNKRDVIGLA QTGSGKTAAF ALPVIQELWN NPSPFFAVVL APTRELAYQI SEQFEAIGGS IGVRSSVIVG GMDMVTQAVA LSKKPHVLVC TPGRLMDHLE NTKGFSLKNL KYLIMDEADR LLDMDFGPII DKILKIIPHE RRTLLFSATM TSKVEKLQRA SLHQPVRVAV SSKFSTVDTL IQRYLFFPFK HKDTYLVYLV NELAGNSIII FARTVNDTQR LAILLRTLGF SAIPLHGQLS QSNRLGALNK FKSGARSTLV ATDVAARGLD IPLVDVVINY DIPTDSKAYI HRVGRTARAG RAGKSIALVT QYDLEPFLRI EATIGKKMQE YEIDKEGVFL LSERVGEAQR EAIQMKEIH DRRKSKGKLH TKRKRDDLDR EEQIY
Source	E.coli
Target Names	rrp3
Protein Names	Recommended name: ATP-dependent rRNA helicase rrp3 EC= 3.6.4.-
Expression Region	1-465
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.