



# Recombinant *Schizosaccharomyces pombe* Transcriptional adapter 2 (ada2)

<b>Product Code</b>	CSB-MP889238SXV
<b>Abbreviation</b>	ada2
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q9P7J7
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	<i>Schizosaccharomyces pombe</i> (strain 972 / ATCC 24843) (Fission yeast)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MPPQKYHCNV CAQDITRSIH IRCVECVDVD LCIPCFTSGA SLGTHHP SHP YRIIETNSYP IFDENWGADE ELLIDACET LGLGNWADIA DYVGNARTKE ECDRHLYKTY IESDCYPLAS VELPGPVDRI AFAARKRARI EAFQPPPIIP QKPLASTPQC HEIQGYMPGR LEFDQEYMNE AELPIKDMNF DDDLHESAKH EMQLKLTMLN IYNSRLTRRA VRKQTIFNHN LLDYRRLQAN EKRM SKEERN LLNKTKAFAR LLTGPDYQKF VNSYHEQITL KKQISDLQEW RQMGLTTLEQ GHKYERDKTQ KFLLSKASAS YDKQLRHVKS FNQTT SAPFQ VRDIQKIVPR KPATPTMFSA SADRQLLED EQALCSKLQI FPKPFLALKF ALISASLTSK KPFQKTDVN LFKHLDANKV EQVYDFFHNA RWIGAPT
<b>Source</b>	Mammalian cell
<b>Target Names</b>	ada2
<b>Protein Names</b>	Recommended name: Transcriptional adapter 2
<b>Expression Region</b>	1-437
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.