



# Recombinant *Schizosaccharomyces pombe* Uncharacterized kinase C24B10.02c (SPCC24B10.02c)

<b>Product Code</b>	CSB-BP889239SXV
<b>Abbreviation</b>	SPCC24B10.02c
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q9P7K3
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	<i>Schizosaccharomyces pombe</i> (strain 972 / ATCC 24843) (Fission yeast)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	MATKVEEINK LKKKNLHFRQ CSADSLSNSS MQFDLSGLNT PVGELATSNL PSPAHPFGE LHQESRTSNS SAMHIENVVA SRLMYNEVAN GSFALPKNI LVVTKPRKHS LVYKTAEITK YILTIGTPET KVVYVDMRLAR SKRFSAHNIA KEANTDIDRI KYWNPYICLI KPSIFDLAIT IGDNSTLLYT SWLFQKIGPP VLSFSDDDVP GFLTHFSLSN YQQHLYQVLT QNVSLRFCSR LQCSFHKYDE KTKQYSLAST TYSLDEILIS RGEHPFISNL NVYNNSELMT VVQADGLVVA TPTGSTNISA NAGGSLVHPA LNAILVTPVC PHTLSFRPII LPDYNVLNVE IPLDSRSSAF FSVDRHESVE MHRGDYLSIV TSHYPFTTIQ NPGYQWTKVL EDKFNWNVRE RQKPF SRKPS LSDVKDTSDD KFDITDNSYC REFSADIDG
<b>Source</b>	Baculovirus
<b>Target Names</b>	SPCC24B10.02c
<b>Protein Names</b>	Recommended name: Uncharacterized kinase C24B10.02c EC= 2.7.1.-
<b>Expression Region</b>	1-449
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.