



Recombinant Schizosaccharomyces pombe Uncharacterized kinase C24B10.02c (SPCC24B10.02c)

Product Code	CSB-EP889239SXV
Abbreviation	SPCC24B10.02c
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9P7K3
Product Type	Recombinant Protein
Immunogen Species	Schizosaccharomyces pombe (strain 972 / ATCC 24843) (Fission yeast)
Purity	≥85% (SDS-PAGE)
Sequence	MATKVEEINK LKKKNLHFRQ CSADSLSNSS MQFDLSGLNT PVGELATSNL PSPAHPFGE LHQESRTSNS SAMHIENVVA SRLMYNEVAN GSFALPKNI LVVTKPRKHS LVYKTAEITK YILTIGTPET KVVYVDMRLAR SKRFSAHNIA KEANTDIDRI KYWNPYICLI KPSIFDLAIT IGDNSTLLYT SWLFQKIGPP VLSFSDDDVP GFLTHFSLSN YQQHLYQVLT QNVSLRFCSR LQCSFHKYDE KTKQYSLAST TYSLDEILIS RGEHPFISNL NVYNNSELMT VVQADGLVVA TPTGSTNISA NAGGSLVHPA LNAILVTPVC PHTLSFRPII LPDYNVLNVE IPLDSRSSAF FSVDRHESVE MHRGDYLSIV TSHYPFTTIQ NPGYQWTKVL EDKFNWNVRE RQKPF SRKPS LSDVKDTSDD KFDITDNSYC REFSADIDG
Source	E.coli
Target Names	SPCC24B10.02c
Protein Names	Recommended name: Uncharacterized kinase C24B10.02c EC= 2.7.1.-
Expression Region	1-449
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.