



Recombinant Schizosaccharomyces pombe VEL1-related protein AC977.05c (SPAC977.05c)

Product Code	CSB-EP889249SXV-B
Abbreviation	SPAC977.05c
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9P7U4
Product Type	Recombinant Protein
Immunogen Species	Schizosaccharomyces pombe (strain 972 / ATCC 24843) (Fission yeast)
Purity	>85% (SDS-PAGE)
Sequence	IRF NLTDLECSRL RGPHCPTYLL KVVGTNATYV GEKSFGLDA LTESKGEFFQ RMLEQEPRLI PRLFTIAEND TANFTPLTFT TYLKTCNPQS IENAMIPFVN TVTSEISFDA WAYTAQNSSR ITGLSNQLMN STLYNVQVAT CTPGFSALLL DSPTINVFNN EEGMPSWCQP IELIPVCPLD EGFN
Source	E.coli
Target Names	SPAC977.05c
Protein Names	Recommended name: VEL1-related protein AC977.05c
Expression Region	18-204
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.