



# Recombinant Mouse Interleukin-17B (IL17b)

<b>Product Code</b>	CSB-EP889559MO-B
<b>Abbreviation</b>	IL17b
<b>Storage</b>	<p>The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.</p> <p>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.</p>
<b>Uniprot No.</b>	Q9QXT6
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Mus musculus (Mouse)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	RNTKGKRK GQGRPSPLAP GPHQVPLDLV SRVKPYARME EYERNLGEMV AQLRNSSEPA KKKCEVNLQL WLSNKRSLSP WGYSINHDPS RIPADLPEAR CLCLGCVNPF TMQEDRSMVS VPVFSQVPVR RRLCPQPPRP GPCRQRVME TIAVGCTCIF
<b>Source</b>	E.coli
<b>Target Names</b>	IL17b
<b>Protein Names</b>	Recommended name: Interleukin-17B Short name= IL-17B Alternative name(s): Cytokine CX1 Cytokine-like protein ZCYTO7 Neuronal interleukin-17-related factor
<b>Expression Region</b>	23-180
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Target Details</b>	This protein is a T cell-derived cytokine that shares sequence similarity with IL17. This cytokine was reported to stimulate the release of TNF alpha (TNF) and IL1 beta (IL1B) from a monocytic cell line. Immunohistochemical analysis of several nerve tissues indicated that this cytokine is primarily localized to neuronal cell bodies.
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	<p>The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.</p> <p>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.</p>