



Recombinant Arabidopsis thaliana F-box protein At1g11810 (At1g11810)

Product Code	CSB-BP890290DOA
Abbreviation	At1g11810
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9SA94
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	≥85% (SDS-PAGE)
Sequence	MTTTMSTLPV VLVDEILARV PITSLRSLRS TCKKWEASSK TNLVGGKATA RKSSHVGFIL IGDKICSMKL DLNGGDDFVD TSVNQVSAFD GFAISQLFHC DGLLFCISNN HYSNYTLMVC NMYLGETRLI QNRSLFESYQ NFCSYAFGYD SSKNRNHKIL RNNSVSGGYE IYSLKSDSWK DLNVDLEKSI HLWRLGSVSL KGNAYFRVIK VIEEGVWEYN LLCFDFTRES FGKLLSLPFE SLGDEGEGIM VISCVRDDHL AVLYQRDTLG IWISTEIEPN KVSWREFLQV DLATLDGFPD VFIAGR FIVD EEKQVVVVFV QETELDLNHG NAFIFGRDGY FTSFTVGDAP PEDFTSYVPS LVSLQIDKTG KRKARDD
Source	Baculovirus
Target Names	At1g11810
Protein Names	Recommended name: F-box protein At1g11810
Expression Region	1-377
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.