



Recombinant Arabidopsis thaliana Zinc finger AN1 and C2H2 domain-containing stress-associated protein 13 (SAP13)

Product Code	CSB-EP890305DOA-B
Abbreviation	SAP13
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9SCM4
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Sequence	MGTPEFPDLG KHCSVDYCKQ IDFLPFTCDR CLQVYCLDHR SYMKHDCPKG NRGDVTVIC PLCAKGVRLN PDEDPNITWE KHVNTDCDPS NYEKAVKKKK CPVPRCRELL TFSNTIKCRD CSIDHCLKHR FGPDHSCSGP KKPESSFSFM GFLSTNTKEA PASSSSSSRW SSLFASAEAS ISRLGNDISQ KLQFASGNDG NSEKTQERNG KQNCGKVTVD VCPKCSRGFR DPVDLLKHID KDHRGTSKA
Source	E.coli
Target Names	SAP13
Protein Names	Recommended name: Zinc finger AN1 and C2H2 domain-containing stress-associated protein 13 Short name= AtSAP13
Expression Region	1-249
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.