



# Recombinant Arabidopsis thaliana Putative U-box domain-containing protein 58 (PUB58)

<b>Product Code</b>	CSB-YP890334DOA
<b>Abbreviation</b>	PUB58
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q9SGT2
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Arabidopsis thaliana (Mouse-ear cress)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	MVENSYVLFA RLCVELFHFL PPLPSDESHG DITTFERLFL RKCRIELENA SPKTEHLPLV YVDETNYRFR IKTVRVLAEV YKNSKITETT RKSMIQVLMN PILPPERITD AMNLFRSIIG KLADFHFSDE KFNQLVRSSR VVELEGNYNE EVKLRKEAED ALAMKKEDVE MMEQLLESYK EEQGKLQLQA KALEHKLEAE LRHRKETETL LAIERDRIEK VKIQLETVEN EIDNTRLKAE EFERKYEGEM ILRRESEIAL EKEKKELEEV KLKLETYERE QENLSSEVRT WQDKYEQESS LRKLSEYALS REQEELQIVK GLLEFYNGEA DAMREERDKA LKTAKEQMEK RQPPSSFFCP ITQEVMKDPH FAADGFTYEA ESIRKWLSTG HQTSPMTNLR LSHLTLV PNR ALRSAIEELV
<b>Source</b>	Yeast
<b>Target Names</b>	PUB58
<b>Protein Names</b>	Recommended name: Putative U-box domain-containing protein 58 EC= 6.3.2.- Alternative name(s): Plant U-box protein 58
<b>Expression Region</b>	1-420
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.