



# Recombinant Arabidopsis thaliana Purple acid phosphatase 10 (PAP10)

<b>Product Code</b>	CSB-EP890353DOA-B
<b>Abbreviation</b>	PAP10
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q9SIV9
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Arabidopsis thaliana (Mouse-ear cress)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	NGGIT SRYVRKLEAT VDMPLDSDVF RVPCGYNAPQ QVHITQGDVE GKAVIVSWVT QEAKGSNKVI YWKENSTKKH KAHGKTNTYK FYNYTSGFIH HCPIRNLEYD TKYYYVLGVG QTERKFWFFT PPEIGPDVPY TFGLIGDLGQ SYDSNITLTH YENNPTKGQA VLFVGDISYA DTYPDHDNRR WDSWGRFAER STAYQPWIWT TGNHELDFAE EIGENRPFKP FTHRYRTPYR SSGSTEPFWY SIKRGPAYII VLASYSAYGK YTPQYQWLEE EFPKVNRTET PWLIVLMHSP WYNSYDYHYM EGETMRVME AWFVKYKVDV VFAGHVHAYE RSERVSNIAY NVVNGICTPV KDQSAPVYIT IGDGGNIEGL ATKMTEPQPK YSAFREASFG HAIFSIKNRT HAHYGWHRNH DGYAVEGDRM WFYNRFWHPV DDSPSCNS
<b>Source</b>	E.coli
<b>Target Names</b>	PAP10
<b>Protein Names</b>	Recommended name: Purple acid phosphatase 10 EC= 3.1.3.2
<b>Expression Region</b>	26-468
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.