



Recombinant Arabidopsis thaliana Putative 60S ribosomal protein L13-2 (RPL13C)

Product Code	CSB-BP890391DOA
Abbreviation	RPL13C
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9SMT4
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Sequence	MKHNNVIPNG HFKKKWENYV KTSFNQPAMK TRRRIARQNK AVKIFPRPTA GPIRPVVHAQ TLTYNMKVRA GKGFTLEELK AAGIPKKLAP TIGISVDHHR KNRSLEGFQT NVQRLKTYKA KLVIFPRCAR TVKVGDSAQQ ELANATQVQV DHMPIVREMP TMELVKLTSD MKLFNAYDKI RLEGINKRHA GVRAKRAAEA EKEEEK
Source	Baculovirus
Target Names	RPL13C
Protein Names	Recommended name: Putative 60S ribosomal protein L13-2
Expression Region	1-206
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.