



Recombinant Arabidopsis thaliana Glutathione S-transferase U8 (GSTU8)

Product Code	CSB-MP890407DOA
Abbreviation	GSTU8
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9SR36
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Sequence	MNQEEHVKLL GLWGSPFSKR VEMVLKCLKGI PYEYIEEDVY GNRSPMLLKY NPIHKKVPVL IHNGRSIAES LVIVEYIEDT WKTTHILPQ DPYERAMARF WAKYVDEKVM LAVKKACWGP ESEREKEVKE AYEGLKCLEK ELGDKLFFGG ETIGFVDIAA DFIGYWLGIF QEASGVTIMT AEEFPKLQRW SEDFVGNNFI KEVLPPKEKL VAVLKAMFGS VTSN
Source	Mammalian cell
Target Names	GSTU8
Protein Names	Recommended name: Glutathione S-transferase U8 Short name= AtGSTU8 EC= 2.5.1.18 Alternative name(s): GST class-tau member 8
Expression Region	1-224
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.