



Recombinant Arabidopsis thaliana Alanine--glyoxylate aminotransferase 2 homolog 3, mitochondrial (At3g08860)

Product Code	CSB-MP890409DOA
Abbreviation	At3g08860
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9SR86
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Sequence	SSVRETETE TKLPKMPPFN YSPPPYDGPS TAEIIAKRRE FLSPALFHFY NTPLNIVEAK MQYVFDENGR RYLDAFGGIA TVSCGHCHPE VVNSVVKQLK LINHSTILYL NHTISDFAEA LVSTLPGDLK VVFFTNSGTE ANELAMMMAR LYTGCDIVS LRNSYHGNA A ATMGATAQSN WKFNVVQSGV HHAINPDPYR GIFGSDGEKY ASDVHDLIQF GTSGQVAGFI GESIQGVGGI VELAPGYLPA AYDIVRKAGG VCIADDEVQSG FARTGTHFWG FQSHGVIPDI VTMAKGIGNG IPLGAVVTPP EIAGVLSRRS YFNTFGGNPM CTAAGHAVLR VLHEEKLQEN ANLVGSHLKR RLTLKKNKYE LIGDVRGRGL MLGVEFVKDR DLKTPAKAET LHLMDQMKEM GVLVGKGGFY GNVFRITPPL CFTLSADDFL VDVMDHAMSK M
Source	Mammalian cell
Target Names	At3g08860
Protein Names	Recommended name: Alanine--glyoxylate aminotransferase 2 homolog 3, mitochondrial EC= 2.6.1.44 Alternative name(s): Beta-alanine-pyruvate aminotransferase 3
Expression Region	32-481
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.