



# Recombinant Arabidopsis thaliana 26S protease regulatory subunit 7 homolog B (RPT1B)

<b>Product Code</b>	CSB-MP890422DOA
<b>Abbreviation</b>	RPT1B
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q9SSB4
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Arabidopsis thaliana (Mouse-ear cress)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MIIRDIMEDP ENEGPHNMRS DFEPLLSVLR LRDYETGDII EFDSTEASEE ALSDITEFGS TEASEAHFEE PYSARIKKVE KEINELAEKI CNLGIKESDT GLAPPNQWDL VSDKQMMQEE QPLL VATCTQ IISPNTEDAK YVVDIKKIGK YVVG LGDKAS PTDIEAGMRV GVDQKKYQIQ IPLPPKIDPS VTMMTVEEKP DATYSDIGGC KEQIEKIREV VELPMLHPEK FVRLGIDPPK GVLCYGPPGS GKTLVARAVA NRTGACFIRV VGSELVQKYI GEGARMVREL FQMARSKKAC ILFFDEIDAI GGARFDDGVG SDNEVQRTML EILYQLDGF ARGNIKVLMA TNRPDILDPA LLRPGRLDRK VEFCLPDLEG RTQIFKIHTR TMSCERDIRF ELLAGLCPNS TGADIRSVCI EAGMYAIGAR RKS VTEKDFL DAVNKVVKGY QKFSATPKYM AYYI
<b>Source</b>	Mammalian cell
<b>Target Names</b>	RPT1B
<b>Protein Names</b>	Recommended name: 26S protease regulatory subunit 7 homolog B Alternative name(s): 26S proteasome AAA-ATPase subunit RPT1b 26S proteasome subunit 7 homolog B Regulatory particle triple-A ATPase subunit 1b
<b>Expression Region</b>	1-464
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.