



Recombinant Arabidopsis thaliana Pentatricopeptide repeat-containing protein At4g32450, mitochondrial (PCMP-H63)

Product Code	CSB-MP890450DOA
Abbreviation	PCMP-H63
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9SUU7
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	≥85% (SDS-PAGE)
Sequence	QNWQSSDGCS SYGTTGNGVP QENNTGGNHF QQDHSGHSSL DELDSICREG KVKKAVEIHK SWRNEGYVVD LPRLFWIAQL CGDAQALQEA KVVHEFITSS VGISDISAYN SIIEMYSGCC SVEDALTVFN SMPERNLETW CGVIRCFANK GQGEDAIDTF SRFKQEGNKP DGEMFKEIFF ACGVLGDMNE GLLHFESMYK EYGIIPCMEH YVSLVKMLAE PGYLDEALRF VESMEPNVDL WETLMNLSRV HGDLILGDRC QDMVEQLDAS RLNKEKAGL VPKSSDLVK EKLQRMAKGP NYGIRYMAAG DISRPENREL YMALKSLKEH MIEIGYVPLS KLALHDVDQE SKDENLFNHN ERFAFISTFL DTPARSLIRV MKNLRVCADC HNALKLMSKI VGRELISRDA KRFHMKDGV CSCREYW
Source	Mammalian cell
Target Names	PCMP-H63
Protein Names	Recommended name: Pentatricopeptide repeat-containing protein At4g32450, mitochondrial
Expression Region	111-537
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients,



storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.