



Recombinant Arabidopsis thaliana Peroxisome biogenesis factor 10 (PEX10)

Product Code	CSB-BP890480DOA
Abbreviation	PEX10
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9SYU4
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Sequence	MRLNGDSGPG QDEPGSSGFH GGIRRFPLAA QPEIMRAAEK DDQYASFIHE ACRDAFRHLF GTRIALAYQK EMKLLGQMLY YVLTGSGQQ TLGEEYCDII QVAGPYGLSP TPARRALFIL YQTAVPYIAE RISTRAATQA VTFDESDEFF GDSHIHSPRM IDLPSSSQVE TSTSVVSRNLN DRLMRSWHRA IQRWPVVLVP AREVLQLVLR ANLMLFYFEG FYYHISKRAS GVRYVFIGKQ LNQRPRYQIL GVFLLIQLCI LAAEGLRRSN LSSITSSIIQ ASIGSYQTSG GRGLPVLNEE GNLITSEAEK GNWSTSDSTS TEAVGKCTLC LSTRQHPTAT PCGHVFCWSC IMEWCNEKQE CPLCRTPNTH SSLVCLYHSD F
Source	Baculovirus
Target Names	PEX10
Protein Names	Recommended name: Peroxisome biogenesis factor 10 Alternative name(s): PER10 Peroxin-10 Peroxisomal biogenesis factor 10 Short name= AtPEX10 Short name= AthPEX10 Peroxisome assembly protein 10 Pex10p
Expression Region	1-381
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.