



Recombinant Drosophila melanogaster J domain-containing protein (jdp)

Product Code	CSB-EP890609DLU
Abbreviation	jdp
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9TVP3
Product Type	Recombinant Protein
Immunogen Species	Drosophila melanogaster (Fruit fly)
Purity	>85% (SDS-PAGE)
Sequence	MSAVDAIINY KRSPNEDFYG LLHCDENSSP EQIQAEYKVL ALQYHPDKNS GDKEAEAKFQ QLKEAKETLC DPEKRAIYDK WRNSGISMSY KQWLGGMKEHV GQVRRYIAE KVDKESMHWV TPKTKDRMLP ETGGAAAQGP GTGAGCSSGG LTAASPNPAH RRASEGGAAL YYGSRKGEWG GENTDVANKF RNYEI
Source	E.coli
Target Names	jdp
Protein Names	Recommended name: J domain-containing protein Alternative name(s): dJDP
Expression Region	1-195
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.