



# Recombinant Drosophila melanogaster Eukaryotic translation initiation factor 3 subunit H (eIF-3p40)

|                          |   |
|--------------------------|---|
| <b>Product Code</b>      | CSB-YP890649DLU   |
| <b>Abbreviation</b>      | eIF-3p40  |
| <b>Storage</b>           | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.  |
| <b>Uniprot No.</b>       | Q9U9Q4  |
| <b>Product Type</b>      | Recombinant Protein   |
| <b>Immunogen Species</b> | Drosophila melanogaster (Fruit fly)   |
| <b>Purity</b>            | ≥85% (SDS-PAGE)   |
| <b>Sequence</b>          | MANRANRHAA RTEDSDNTIN YVQCDGLAVM KMKVHCHEES<br>SNMDLAQGAL LGLVVDKCLE ITNCFPPPKS GDETMDEEMY QLTVMRRLLR<br>VNVDHLHVGW YQSSDVGNSL SMALLESQYH YQTSIEESV VVYDTQKSSR<br>GFLCLKAYRL TPQAIQMYKD GDFTPEAFRT LKVGYENLFA EIPIVIKNSP<br>LTNIMMSELN ELLPEDKGFHN FLDLGTATVL ENQMRS LIER VDELYQEAVR<br>YNKYQQVVK QDTEKHRALA KLAENAVRT SKGEPTVPEE EVIKQFRPMT<br>APNRLTATIT SGQINTHAQH IAQFCSQSLA KLFITESLQI AKEAKETK        |
| <b>Source</b>            | Yeast   |
| <b>Target Names</b>      | eIF3h   |
| <b>Protein Names</b>     | Recommended name: Eukaryotic translation initiation factor 3 subunit H Short name= eIF3h Alternative name(s): Eukaryotic translation initiation factor 3 subunit 3  |
| <b>Expression Region</b> | 1-338   |
| <b>Notes</b>             | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.   |
| <b>Tag Info</b>          | Tag type will be determined during the manufacturing process.   |
| <b>Protein Length</b>    | full length protein   |
| <b>Reconstitution</b>    | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference. |
| <b>Shelf Life</b>        | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.  |