



Recombinant Human Serine/threonine-protein kinase NLK (NLK)

Product Code	CSB-MP890657HU
Abbreviation	NLK
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9UBE8
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MSLCGARANA KMMAAYNGGT SAAAAGHHHH HHHHLPHLPP PHLHHHHHPQ HHLHPGSAAA VHPVQQHTSS AAAAAAAAAA AAAMLNPGQQ QPYFSPAPG QAPGPAAP AQQVAAAAAT VKAHHHQHSH HPQQQLDIEP DRPIGYGAFG VVWSVTDPRD GKRVALKKMP NVFQNLVSCK RVFRELKMLC FFKHDNVLSA LDILQPPHID YFEEIYVVTE LMQSDLHKII VSPQPLSSDH VKVFLYQILR GLKYLHSAGI LHRDIKPGNL LVNSNCVLKI CDFGLARVEE LDESRHMTQE VVTQYYRAPE ILMGSRHYSN AIDIWSVGC FAELLGRRIL FQAQSPIQQL DLITDLLGTP SLEAMRTACE GAKAHILRGP HKQPSLPVLY TLSSQATHEA VHLLCRMLVF DPSKRISAKD ALAHPYLDEG RLRYHTCMCK CCFSTSTGRV YTSDFEPVTN PKFDDTFEKN LSSVRQVKEI IHQFILEQQK GNRVPLCINP QSAAFKSFIS STVAQPSEMP PSPLVWE
Source	Mammalian cell
Target Names	NLK
Protein Names	Recommended name: Serine/threonine-protein kinase NLK EC= 2.7.11.24 Alternative name(s): Nemo-like kinase Protein LAK1
Expression Region	1-527
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at $-20^{\circ}\text{C}/-80^{\circ}\text{C}$. The shelf life of lyophilized form is 12 months at $-20^{\circ}\text{C}/-80^{\circ}\text{C}$.