



Recombinant Human SRR1-like protein (SRRD)

Product Code	CSB-EP890683HU-B
Abbreviation	SRRD
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9UH36
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MAAAAAAAL SWQAAAPRKR RSAARRPRRR EAAPRGREAA PRGREAAPRG PEAEFESDSG VVLRRIWEAE KDLFISDFWS SALETINRCL TKHLEQLKAP VGTLSDFGN LHLDSLPEES DVATDSIPRE ILVTGTCHLK CVCYGIGNFA TCIVARNQLT FLLLLLEKCQ IPRSHCWVYD PLFSQLEIEV LNTLGVTVLS ENEEGKRSIR GEPTIFYMLH CGTALYNNLL WSNWSVDALS KMVIIGNSFK GLEERLLARI LQKNYPYIAK ILKGLEEELF PQTSQYMDIF NDTSVHWFPV QKLEQLSIDI WEFREEPDYQ DCEDLEIIRN KREDPSATD
Source	E.coli
Target Names	SRRD
Protein Names	Recommended name: SRR1-like protein Alternative name(s): SRR1 domain-containing protein
Expression Region	1-339
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.