



Recombinant Human Regulating synaptic membrane exocytosis protein 3 (RIMS3)

Product Code	CSB-EP890726HU-B
Abbreviation	RIMS3
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9UJD0
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MFNGEPGPAS SGASRNVVRS SSISGEICGS QQAGGGAGTT TAKKRRSSLG AKMVAIVGLT QWSKSTLQLP QPEGATKKLR SNIRRSTETG IAVEMRSRVT RQGSRESTDG STNSNSSDGT FIFPTTRLGA ESQFSDFLDG LGPAQIVGRQ TLATPPMGDV HIAIMDRSGQ LEVEVIEARG LTPKPGSKSL PATYIKVYLL ENGACLAKKK TKMTKKKCDP LYQQALLFDE GPQGKVLQVI VWGDYGRMDH KCFMGMAQIM LDELDSA AV TGWYKLFPTS SVADSTLGSL TRRLSQSSLE SATSPSCS
Source	E.coli
Target Names	RIMS3
Protein Names	Recommended name: Regulating synaptic membrane exocytosis protein 3 Short name= Nim3 Alternative name(s): RIM3 gamma Rab-3-interacting molecule 3 Short name= RIM 3
Expression Region	1-308
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.