



Recombinant Human Cleavage and polyadenylation specificity factor subunit 3 (CPSF3)

Product Code	CSB-BP890734HU
Abbreviation	CPSF3
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9UKF6
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	SAIPAEESD QLLIRPLGAG QEVGRSCIIL EFKGRKIMLD CGIHPGLEGM DALPYIDLID PAEIDLLIS HFHLDHCGAL PWFLQKTSFK GRTFMTHATK AIYRWLLSDY VKVSNISADD MLYTETDLEE SMDKIETINF HEVKEVAGIK FWCYHAGHVL GAAMFMIEIA GVKLLYTADF SRQEDRHLMA AEIPNIKPDI LIIESTYGTH IHEKREEREA RFCNTVHDIV NRGGRGLIPV FALGRAQELL LILDEYWQNH PELHDIPIYY ASSLAKKCMA VYQTYVNAMN DKIRKQININ NPFVFKHISN LKSMDFHDDI GPSVVMASPG MMQSGLSREL FESWCTDKRN GVIIAGYCVE GTLAKHIMSE PEEITMSGQ KLPLKMSVDY ISFSAHTDYQ QTSEFIRALK PPHVILVHGE QNEMARLKA LIREYEDNDE VHIEVHNPRN TEAVTLNFRG EKLAKVMGFL ADKKPEQGQR VSGILVKRNF NYHILSPCDL SNYTDLAMST VKQTQAIPYT GPFNLLCYQL QKLTGDVEEL EIQEKPAKV FKNITVIQEP GMVVLEWLAN PSNDMYADTV TTVILEVQSN PKIRKGAVQK VSKKLEMHVY SKRLEIMLQD IFGEDCVSVK DDSILSVTV D GKTANLNLET RTVECEEGSE DDESLREMVE LAAQRLYEAL TPVH
Source	Baculovirus
Target Names	CPSF3
Protein Names	Recommended name: Cleavage and polyadenylation specificity factor subunit 3 EC= 3.1.27.- Alternative name(s): Cleavage and polyadenylation specificity factor 73 kDa subunit Short name= CPSF 73 kDa subunit mRNA 3'-end- processing e
Expression Region	2-684
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein



Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.