



Recombinant Human TBC1 domain family member 24 (TBC1D24)

Product Code	CSB-BP890752HU
Abbreviation	TBC1D24
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9ULP9
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MDSPGYNCFV DKDKMDAAIQ DLGPKELSCT ELQELKQLAR QGYWAQSHAL RGKVYQRLIR DIPCRTVTPD ASVYSDIVGK IVGKHSSSCL PLPEFVDNTQ VPSYCLNARG EGAVRKILLC LANQFPDISF CPALPAVVAL LLHYSIDEAE CFEKACRILA CNDPGRRLLID QSFLAFESSC MTFGDLVNKY CQAAHKLMVA VSEDVLQVYA DWQRWLFGEPLCYFARVFD VFLVEGYKVL YRVALAILKF FHKVRAGQPL ESDSVKQDIR TFVRDIAKTV SPEKLEKAF AIRLFSRKEI QLLQMANEKA LKQKGITVKQ KSVSLSKRQF VHLAVHAENF RSEIVSVREM RDIWSWVPER FALCQPLLLF SSLQHGYSLA RFYFQCEGHE PTLILLIKTTQ KEVCGAYLST DWSERNKFGG KLGFFGTGEC FVFRQLQPEVQ RYEWVVIKHP ELTKPPPLMA AEPTAPLSHS ASSDPADRLS PFLAARHFNL PSKTESMFMA GGSDCLIVGG GGGQALYIDG DLNRGRTSHC DTFNNQPLCS ENFLIAAVEA WGFQDPDTQ
Source	Baculovirus
Target Names	TBC1D24
Protein Names	Recommended name: TBC1 domain family member 24
Expression Region	1-559
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.