



Recombinant Human Forkhead box protein J3 (FOXJ3)

Product Code	CSB-MP890778HU
Abbreviation	FOXJ3
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9UPW0
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MGLYGQACPS VTSLRMTSEL ESSLTSMDWL PQLTMRAAIQ KSDATQNAHG TGISKKNALL DPNTTLDQEE VQQHKDGKPP YSYASLITFA INSSPKKKMT LSEIQWICD NFPYYREAGS GWKNSIRHNL SLNKCFLKVP RSKDDPGKGS YWAIDTNPKE DVLPTRPKKR ARSVERASTP YSIDSDSLGM ECIISGSASP TLAINVTNK VTLYNTDQDG SDSPRSSLNN SLSDQSLASV NLNSVGSVHS YTPVTSHPES VSQSLTPQQQ PQYNLPERDK QLLFSEYNFE DLSASFRSLY KSVFEQSLSQ QGLMNIPSES SQQSHTSCTY QHSPSSTVST HPHSNQSSLS NSHGSGLNNTT GSNSVAQVSL SHPQMHTQPS PHPPHRPHGL PQHPQRSPPH APHPQQHSQL QSPHPQHPSH HQHIQHHPNH QHQTTLTHQAP PPPQQVSCNS GVSNDWYATL DMLKESCRIS SSVNWSVDVL SQFQGLMESM RQADLKNWSL DQVQFADLCS SLNQFFTQTG LIHSQSNVQQ NVCHGAMHPT KPSQHIGTGN LYIDSRQNLP PSVMPPPGYP HIPQALSTPG TTMAGHHRAM NQQHMMPSQA FQMRRSLPPD DIQDDFDWDS IV
Source	Mammalian cell
Target Names	FOXJ3
Protein Names	Recommended name: Forkhead box protein J3
Expression Region	1-622
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.