



Recombinant Schizosaccharomyces pombe SH3 domain-containing protein PJ696.02 (SPAPJ696.02)

Product Code	CSB-MP890786SXV
Abbreviation	SPAPJ696.02
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9URW6
Product Type	Recombinant Protein
Immunogen Species	Schizosaccharomyces pombe (strain 972 / ATCC 24843) (Fission yeast)
Purity	>85% (SDS-PAGE)
Sequence	MGLHNPLPSS LKSECKKAGK ILTSFVDPRQ TLGAQEVIPP SVLTNAKGLV IMTVLKAGFL FSGRIGSGLI VARLDDGTWS APSAVMTGGM GVGAQIGSEL TDFVIILNSK AAVQTFARLG SITLGGNLSI AAGPLGRNAE AGGGASVGGM APMFSYSKTK GLFAGVSLEG SVLVERRDAN RSLYRGDITA KRLLSGQVAQ PAAADPLYRV LNSKIFNLNR GDEGDIYNDV PIYADDEPED IWGPSSKSTK RRDSADRSSS YSRRGDSYRS NRSRAHDDDD EDDYSFSRSK SLSRKTAGGS LRSSKMDNRR SKYADTPSPR RSRYSDEDE ESVYSSDVST ESSSQFSSRS SEYSKPSRPT APKPKFKQDS LGPNQARAMY SFAGEQPGDL SFQKGDII DI VERSGSHDDW WTGRIGYREG IFPANYVKLS
Source	Mammalian cell
Target Names	SPAPJ696.02
Protein Names	Recommended name: SH3 domain-containing protein PJ696.02
Expression Region	1-430
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.