



Recombinant Schizosaccharomyces pombe CTD kinase subunit gamma (ctk3)

Product Code	CSB-EP890797SXV
Abbreviation	ctk3
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9USJ8
Product Type	Recombinant Protein
Immunogen Species	Schizosaccharomyces pombe (strain 972 / ATCC 24843) (Fission yeast)
Purity	>85% (SDS-PAGE)
Sequence	MDPFEGRMTF LQLLGKLNAS QFSQIKPAQF AIKHLBLEED LYSCIWEELE SGSFNTRVNI MYFVDTLCEM CLKNGLTGGY LNMISRDICK LVQNVAPIGA AGAANAPEVR KVLQSLHEKK VIDDNQYKDA MATVEAHEQA SKSGDTSTSG AISKNDILKR IEEDRERHKR MRENIWAISE PELEAEIAWN TTQGITESDL ESLKDEYEKF NECLHATS
Source	E.coli
Target Names	ctk3
Protein Names	Recommended name: CTD kinase subunit gamma Short name= CTDK-I gamma subunit Alternative name(s): CTD kinase subunit 3
Expression Region	1-218
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.