



# Recombinant *Schizosaccharomyces pombe* UPF0652 protein C16A11.03c (SPCC16A11.03c)

<b>Product Code</b>	CSB-MP890801SXV
<b>Abbreviation</b>	SPCC16A11.03c
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q9USN2
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	<i>Schizosaccharomyces pombe</i> (strain 972 / ATCC 24843) (Fission yeast)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MSIAQLFNQS WVEKSKYVPV RLSDEERKIL SLLEAALEVI DYTGHVDIIS YTRTRKLMVK YLNEMCSIVM GLVTAMDIKA GRDLLIGRDH KSNARFFQTV FEIGRRYKIM NPEKMRATYG AVMYMCQDSL IPDVRSQLGF DFVSPIKTVY NVLEKHKLLG LLEEKQLLNN LLKQSDPQSN ANAKAELLKE REEASETLLK KYNPGKDEKL QKVLTD CFAS LADHEAFLLA NRNPVEKMRA YLHKFFNPHD TKNGSLKIGY MTGAKLNHDH KTQFFYVDQS LVFWSCMMDQ MFLLWLES DA SLLDKHSRYF ISDTGQGLNR VQLCPLVRST VTRILSSVQK QQEIPWMGSS VIHLGDRDVP NALMFIDKYR QVPHILAPLV KVLQQLEFLR DPYLVQYIEN EYGSVNGLQK TILLDFFRHG FNGQGS DGGG CIDGRLTSAW NWTNEITKKK YYRILLMSGF LNFEGI
<b>Source</b>	Mammalian cell
<b>Target Names</b>	SPCC16A11.03c
<b>Protein Names</b>	Recommended name: UPF0652 protein C16A11.03c
<b>Expression Region</b>	1-466
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.