



Recombinant Schizosaccharomyces pombe 20S-pre-rRNA D-site endonuclease nob1 (nob1)

Product Code	CSB-EP890829SXV-B
Abbreviation	nob1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9UTK0
Product Type	Recombinant Protein
Immunogen Species	Schizosaccharomyces pombe (strain 972 / ATCC 24843) (Fission yeast)
Purity	>85% (SDS-PAGE)
Sequence	MSKSITHLVL DTGGIICSST LLRNSAESFY TIPRVIAEIR DETSRKNFEL WGDQVIQRVP KPEFIKKVSE FAKQTGDYSS LSVTDIQILA LTYELEVEFN GGDWRLRKYP GQKHINGKPP SNSNSTEDAS KPTSSDTASV KETENSDBPKS AENEVLEGET TQHSNNKEAH PNTEENKEQE DNEEDEDDEG WITPSNIRKK KAEDGVGESL VQPKHLKVAC ATDFSMQNV LLQIGLNLVS SDGFKIQNVK RFVLRCHGCY TVVKDMEKKF CPSCGGNTLI KTTCSINSKG EFQVHLKKNF EWKTRGTKYS LPKPVHGTSN GKGKKNPVLR EDQPEYQRAV RRMQRKKEID LMDEDYLP SL LGVTKDRMY VQIGAGRKNP NEVRRKKR
Source	E.coli
Target Names	nob1
Protein Names	Recommended name: 20S-pre-rRNA D-site endonuclease nob1 EC= 3.1.-.- Alternative name(s): Pre-rRNA-processing endonuclease nob1
Expression Region	1-388
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.