



Recombinant Schizosaccharomyces pombe Uncharacterized peptidase C18A7.01 (SPBC18A7.01, SPBC4F6.19c)

Product Code	CSB-BP890846SXV
Abbreviation	SPBC18A7.01, SPBC4F6.19c
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9UUD8
Product Type	Recombinant Protein
Immunogen Species	Schizosaccharomyces pombe (strain 972 / ATCC 24843) (Fission yeast)
Purity	>85% (SDS-PAGE)
Sequence	MVSFESSFER GTDFLNRNFK KCLFACISIF IFALLALSFL SLLQPDTVQR LYQCAVPSMI YVPPMINEAI SIQHEEFNNR RRRLSAALRE DKLDALIMEP TVSMDYFANI TTGSWGLSER PFLGIIFSDD EPYPGDVASR IYFLVPKFEL PRAKELVGKN IDAKYITWDE DENPYQVLYD RLGPLKLMID GTVRNFIAQG LQYAGFTTFG VSPRVASLRE IKSPAEDDIM SRVNIATVAA IRSVQPCIKP GITEKELAEV INMLFVYGGL PVQESPIVLF GERAAMPHGG PSNRRLLKSE FVLMDVGTTL FGYHSDCTRT VLPHGQKMTE RMEKLWNLVY DAQTAGIQML SHLSNTSCAE VDLAARKVIK DAGYGEYFIH RLGHGLGLEE HEQTYLNPAN KGTPVQKGNV FTVEPGIYIP DEIGIRIEDA VLASDVPILL TNFRAKSPYE P
Source	Baculovirus
Target Names	SPBC18A7.01
Protein Names	Recommended name: Uncharacterized peptidase C18A7.01 EC= 3.4.-.-
Expression Region	1-451
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.