



Recombinant *Schizosaccharomyces pombe* Peptidyl-prolyl cis-trans isomerase cyp6 (cyp6)

Product Code	CSB-MP890847SXV
Abbreviation	cyp6
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9UUE4
Product Type	Recombinant Protein
Immunogen Species	<i>Schizosaccharomyces pombe</i> (strain 972 / ATCC 24843) (Fission yeast)
Purity	≥85% (SDS-PAGE)
Sequence	MSVLIETTVG DLVIDLFFVKE APKTCENFLK LCKLKYYNFC PFYNIQHNYT CQTGDPLGPT GDGGRCVWNV LNKGTRFFKA EFNPSLVHMK MGLVSMSTAT ISSRDDKLLV CGSQFIITLS DNLEGLDERY PIYGQVAEGF DTLLKINDAI CDEEGQPYRD IRIKHTIILD DPFEDPPDLV EPLRSPSPTP EQLATVRIGE NEQIESETSE DKLQREKEME AEAEAVTLEM IGDLPFAHVA PPENLVFVCK LNPVTQDEDL ELIFSRFGKI ISCQVIRDKE TGDSLQYAFI EFDNKESVEK AYFKMQNVLV DDSRIHVDFS QSVARYRQYY NSNRDRKRSS SRSDDREYHR RSDGRYDRSN YRDDYRHRRK ERDHRDDQSS FRNERFSNYY GDDRSYHKRR NTGNKNCDH LRDKSPERRY RYDRRYRDDR YR
Source	Mammalian cell
Target Names	cyp6
Protein Names	Recommended name: Peptidyl-prolyl cis-trans isomerase cyp6 Short name= PPlase cyp6 EC= 5.2.1.8 Alternative name(s): Rotamase cyp6
Expression Region	1-432
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.