



Recombinant Schizosaccharomyces pombe Ribosome production factor 2 homolog (SPAC926.08c)

Product Code	CSB-BP890849SXV
Abbreviation	SPAC926.08c
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9UUG1
Product Type	Recombinant Protein
Immunogen Species	Schizosaccharomyces pombe (strain 972 / ATCC 24843) (Fission yeast)
Purity	>85% (SDS-PAGE)
Sequence	MLRQVKPKNA RTKRALEKRE PKLVEGAKTA IFLRGNATSQ ISLDVLGDIH ALKKPFVSVNF QKKNNILPFE DASSLEFFSE KNDAAALAVMA THNKKRPHNL TWVRRFFNYRV LDMIELGIVN YKSIQSFSA PIVPGTKPMI LFQGPVFDAAH PTYRHIKSLF LDFFRGEPIQ KLDSAGLSYV IVVSAEAQAE DETKPLPLVH FRVYGTKLLK TKTNLPRVEL EEMGPRIDFN IRRVQPAESD VLEEALKKPK TQEPKPKKNV DVDIIGNKVG RIHVDQQDLG NLQTRKMKGL KRSVEEREDS ENEEVEIEED VISDASE
Source	Baculovirus
Target Names	SPAC926.08c
Protein Names	Recommended name: Ribosome production factor 2 homolog Alternative name(s): Brix domain-containing protein 1 homolog Ribosome biogenesis protein RPF2 homolog
Expression Region	1-317
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.