



# Recombinant Human Lariat debranching enzyme (DBR1)

<b>Product Code</b>	CSB-EP890914HU-B
<b>Abbreviation</b>	DBR1
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q9UK59
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	MRVAVAGCCH GELDKIYETL ALAERRGPGP VDLLLLCCGDF QAVRNEADLR CMAVPPKYRH MQTFYRYYSG EKKAPVLTFL IGGNHEASNH LQELPYGGWV APNIYYLGLA GVKYRGVRI GGISGIFKSH DYRKGHFECF PYNSSTIRSI YHVRNIEVYK LKQLKQPIDI FLSDHWPRSI YHYGNKKQLL KTKSFFRQEV ENNTLGSPAA SELLEHLKPT YWFS AHLVK FAALMQHQAK DKGQTARATK FLALDKCLPH RDFLQILEIE HDPSAPDYLE YDIEWLTILR ATDDLINVTG RLWNMPENNG LHARWDYSAT EEGMKEVLEK LNHDLKVPCN FSVTAACYDP SKPQTQMQLI HRINPQTTEF CAQLGIIDIN VRLQKSKEEH HVCGEYEEQD DVESNDSGED QSEYNTD TSA LSSINPDEIM LDEEEDSDSI VSAHSGMNTP SVEPSDQASE FSASFSDVRI LPGSMIVSSD DTV DSTIDRE GKPGGTVESG NGEDLTKVPL KRLSDEHEPE QRKKIKRRNQ AIYAAVDDDD DDAA
<b>Source</b>	E.coli
<b>Target Names</b>	DBR1
<b>Protein Names</b>	Recommended name: Lariat debranching enzyme EC= 3.1.-.-
<b>Expression Region</b>	1-544
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.